



SEAWORK MARINE SERVICES LTD

Unit 6, Budds Lane, Romsey, Hampshire, SO51 0HA
 Telephone: +44 (0)8455 192833
 Email: info@seaworkmarineservices.com
www.seaworkmarineservices.com

Seawork are one of a few Inland Inshore contractors that own and operate Decompression Chambers. We currently have four, one is a Containerised Twin Lock Chamber, two are Containerised Twin Lock Chamber with 3 Diver Panels and one is a Twin Lock Chamber with 2 Diver Panel.

This unit (Con 06) is the Containerised Chamber System consisting of Twin Lock Chamber Container with LP and HP air supplies.



This system is housed in a 20 foot ISO Shipping Container that is insulated, lined, heated and has primary and emergency lighting. This unit is equipped with a side personnel access door and window as well as penetrator plates for connection of air supplies and power.

The chamber is a 2 compartment (twin lock) unit with Medlock and is built to BS5500 and is fully tested and approved by Lloyds Register who also witness the 5 yearly hydro test regimes.

Audited to IMCA DO23 with full certification packages and equipped with First Aid kit, DMAC015 medical kit, Defib and an oxygen therapy kit.



Containerised Decompression Chamber – Con 06





SEAWORK MARINE SERVICES LTD

Unit 6, Budds Lane, Romsey, Hampshire, SO51 0HA
Telephone: +44 (0)8455 192833
Email: info@seaworkmarineservices.com
www.seaworkmarineservices.com


Kept in meticulous condition, as you would expect for a Diving Chamber and HBOT Unit by our in house fully qualified Diving Technicians with decades of experience this unit is ready to be shipped to site at immediate notice.

Inside the chamber it is equipped with Primary communications, back up sound powered phones, hyperbaric fire extinguishers, BIBS, Internal Lighting in both locks together with fire retardant bedding, caisson gauge in the inner lock.

Externally there are depth Gauges and Therapeutic gauges, viewports into both locks as well as oxygen and CO2 monitoring with crystal clear communications all housed in a temperature stabilised secure safe environment.

To operate the chamber to its full potential there is a 500 ltr, LP receiver to provide main air to the chamber and a 6 x 65ltr cylinder, 350 bar bank to provide backup reserve air to the chamber all housed and rigidly plumbed into a fully operational unit.

This Containerised Chamber is equipped with 20 metre interconnecting deck hoses for both the main (LP) and the Secondary (HP) air as well as the power supply cable.



Containerised Decompression Chamber – Con 06

